

2002

5' Fluorogenic PCR Assay for the iap Gene of *Listeria* species

Irene V. Wesley

United States Department of Agriculture

Alissa D. Jourdan

United States Department of Agriculture

H. Scott Hurd

United States Department of Agriculture

Angelia L. Pinkerton

United States Department of Agriculture

Follow this and additional works at: http://lib.dr.iastate.edu/swinereports_2001



Part of the [Agriculture Commons](#), and the [Animal Sciences Commons](#)

Extension Number: ASL R1790

Recommended Citation

Wesley, Irene V.; Jourdan, Alissa D.; Hurd, H. Scott; and Pinkerton, Angelia L., "5' Fluorogenic PCR Assay for the iap Gene of *Listeria* species" (2002). *Swine Research Report, 2001*. 26.

http://lib.dr.iastate.edu/swinereports_2001/26

This report is brought to you for free and open access by the Animal Science Research Reports at Iowa State University Digital Repository. It has been accepted for inclusion in Swine Research Report, 2001 by an authorized administrator of Iowa State University Digital Repository. For more information, please contact digirep@iastate.edu.

5' Fluorogenic PCR Assay for the iap Gene of Listeria species

Abstract

We developed a fluorogenic 5' nuclease polymerase chain reaction (PCR) assay targeting the iap gene, which is highly conserved in all Listeria species. We used this assay to screen tissues obtained during a survey of cull sows. To date, we have detected Listeria species in caecal contents (66%, 23 of 35 samples), tonsil homogenates (67%, 4 of 6 samples), and ileocaecal lymph nodes (75%, 6 of 8 samples), but infrequently in feces collected before slaughter (4.4%, 1 of 23 samples). Taken together, these data indicate that second-generation 5' nuclease PCR (TaqMan) assays accelerates detection and quantitation of Listeria during pork production. This may reduce product recalls as well as reduce human illness and deaths. The development of portable instruments for 5' nuclease assays will expedite field surveys.

Keywords

ASL R1790

Disciplines

Agriculture | Animal Sciences

5' Fluorogenic PCR Assay for the *iap* Gene of *Listeria* species.

Irene V. Wesley, Alyssa D. Jourdan, H. Scott Hurd, and Angelia L. Pinkerton, National Animal Disease Center, Ames, IA

ASL-R1790

Summary

We developed a fluorogenic 5' nuclease polymerase chain reaction (PCR) assay targeting the *iap* gene, which is highly conserved in all *Listeria* species. We used this assay to screen tissues obtained during a survey of cull sows. To date, we have detected *Listeria* species in caecal contents (66%, 23 of 35 samples), tonsil homogenates (67%, 4 of 6 samples), and ileocaecal lymph nodes (75%, 6 of 8 samples), but infrequently in feces collected before slaughter (4.4%, 1 of 23 samples). Taken together, these data indicate that second-generation 5' nuclease PCR (TaqMan) assays accelerates detection and quantitation of *Listeria* during pork production. This may reduce product recalls as well as reduce human illness and deaths. The development of portable instruments for 5' nuclease assays will expedite field surveys.

Introduction

Listeria is a nonspore forming, gram-positive, facultatively anaerobic bacillus, that is found in soil, water, vegetation, fish, birds, insects, and a variety of mammals (9). There are six species of *Listeria*: *L. monocytogenes*, *L. ivanovii*, *L. seelegeri*, *L. murrayi*, *L. welshimeri*, and *L. innocua* (9).

L. monocytogenes affects immunocompromised human hosts, including infants and neonates, the elderly, and pregnant women, resulting in meningitis, bacteremia, mild gastroenteritis, and spontaneous abortions (4,9). In the United States, ~2,500 listeriosis cases are reported annually at a medical cost of ~\$0.25 billion. The case fatality rate of 30% is the highest of any foodborne bacterial agent (9). There is "zero tolerance" for *L. monocytogenes* in ready to eat products. From 1990 to 1997 10 of 12 class I (voluntary) recalls involved pork (4). In 1999, 50% of food recalls were due to *L. monocytogenes* contamination (4). *L. innocua*, is usually found in concert with *L. monocytogenes* and thus may be an indicator of its presence (4).

In France, one of the largest listeriosis outbreaks was traced to pickled pork tongue and involved 279 human cases (33% pregnancy related; 7). A second outbreak in France involved 38 patients (82% pregnancy related) and incriminated contaminated pork rillettes (pate). A third episode occurred in spring 2000. Healthy carrier hogs were thought to be the most probable source of *L. monocytogenes* contamination in the plant (5). Because of

its public health significance, highly specific polymerase chain reaction (PCR)-based strategies to identify *L. monocytogenes* have been detailed (1).

The invasion-associated protein (*iap*) gene is conserved in all species of *Listeria* and encodes the p60 protein for the invasion of nonprofessional phagocytic cell (2). Thus, the *iap* gene is an ideal candidate for the detection of all *Listeria* species (6). Herein we report a fluorogenic 5' nuclease PCR assay specific for the *iap* locus.

Materials and Methods

5' Nuclease Assay. Sequences and concentrations for the primers and probe targeting the *iap* gene in the 5' nuclease fluorogenic PCR assay are shown in Table 1.

Bacterial Strain. Strains used to determine specificity were obtained from the National Animal Disease Center (Table 2).

Genomic DNA Extraction. A loop of bacterial growth was harvested from the surface of BHI plates into 200 μ l of sterile water and boiled (15 minutes in a heat block at 103.0°C). The boiled suspensions were pelleted (11,500 rpm) and the supernates were stored (-20°C) for PCR analysis.

Fluorogenic PCR. Each amplification reaction (50 μ l) consisted of template DNA (5 μ l), 5 mM MgCl₂, 0.4 mM dNTPs, 1X AmpliTaq Gold Buffer + Rox, and 2.5 U AmpliTaq Gold DNA polymerase. Amplification conditions included initial denaturation (95°C for 10 minutes) and 35 cycles of denaturation (95°C for 15 seconds), primer annealing (55°C for 30 seconds), and primer extension (72°C for 45 seconds) followed by a final extension (72°C for 3 minutes) and an indefinite hold (25°C). PCR reactions were performed in the ABI Prism 7700 Sequence Detector System and analyzed using the Sequence Detector Software, version 1.6.3 (Applied Biosystems, Foster City, CA).

Cull Sow Tissues: The optimized assay was used to detect *Listeria* spp in tissues collected during a survey of cull sows (8). Fecal samples and tissues were collected from cull sows (n = 181) during five sampling points (January –March 2001). For this study, 72 samples, including feces, caecal contents, ileocaecal lymph nodes, and tonsil homogenates, were screened for *Listeria*.

Results

Fluorogenic 5' Nuclease PCR Assay. The optimized 5' nuclease assay amplified *L. monocytogenes*, *L. innocua*, *L. seelegeri*, *L. ivanovi*, and *L. welshimeri* (Table 2). No amplification was seen with *L. murrayi*.

Genomic DNA extraction. Initially, 5' nuclease assays were conducted using bacteria harvested from Palcam. However, because Palcam contains acriflavine, which interferes with the TaqMan assay, it was necessary to subculture to BHI and then proceed with the DNA extraction.

Cull Sow Tissues. Preliminary analysis of the samples screened to date indicate *Listeria* in caecal contents (66%, 23 of 35 samples), tonsil homogenates (67%, 4 of 6 samples), and ileocaecal lymph nodes (75%, 6 of 8 samples), but infrequently in feces collected before slaughter (4.4%, 1 of 23 samples).

Discussion

The *iap* gene, which encodes the P60 protein needed for penetration into nonphagocytic cells, exists in all species of *Listeria*, including *L. innocua*.

We have designed a TaqMan 5' nuclease assay, which targets the *iap* gene, which is conserved in all species of *Listeria* of public health significance. We used this assay to screen tissues obtained during a survey of cull sows. To date, we have detected *Listeria* species in caecal contents (66%, 23 of 35 samples), tonsil homogenates (67%, 4 of 6 samples), and ileocaecal lymph nodes (75%, 6 of 8 samples), but infrequently in feces collected before slaughter (4.4%, 1 of 23 samples).

Taken together, these data indicate that second generation 5' nuclease PCR (TaqMan) assays accelerate detection and quantitation of *Listeria* during pork production. This may reduce product recalls as well as reduce human illness and deaths. The development of portable instruments for 5' nuclease assays will expedite field surveys.

Acknowledgements

This study was partially funded by the National Alliance for Food Safety and the National Pork Producers Council.

References

1. Bassler, Heidi A., Susan J. Flood, Kenneth J. Livak, Jeffrey Marmaro, Ruth Knorr, and Carl A. Batt. 1995. Use of a fluorogenic probe in a PCR-based assay for the detection of *Listeria monocytogenes*. *Appl. Environ. Microbiol.* 61: 3724 - 3728.
2. Bubert, Andreas, Stephan Kohler, and Werner Goebel. 1992. The homologous and heterologous regions within the *iap* gene allow genus and species specific identification of *Listeria* spp. by polymerase chain reaction. *Appl. Environ. Microbiol.* 58: 2625 - 2632.
3. Buncic S. 1991. The incidence of *Listeria monocytogenes* in slaughtered animals, in meat and in meat products in Yugoslavia. *Int. J. Food Microbiol.* 112: 173 - 180.
4. Farber JM et al. 1998. Incidence and behavior of *Listeria monocytogenes* in meat products. In *Listeria, Listeriosis and Food Safety*, ET Ryser and EH Marth (editors), pp. 505 - 564. Marcel Dekker, Inc., New York.
5. Goulet V. et al. 1993. Listeriosis outbreak associated with consumption of rillettes in France in 1993. *J. Infect. Dis.* 177: 155 - 160.
6. Hein, Ingeborg, Dieter Klein, Angelika Lehner, Andreas Bubert, Ernst Brandl, and Martin Wagner. 2001. Detection and quantification of the *iap* gene of *Listeria monocytogenes* and *Listeria innocua* by a new real-time quantitative assay. *Res. Microbiol.* 152: 37 - 46.
7. J.D. McKean, H.S.Hurd, S. Larsen, R. Griffith, I.Wesley 2001. Impact of commercial pre-harvest processes on the prevalence of salmonella in cull sows. 4th International Symposium on Epidemiology and Control of Salmonella and other Foodborne Pathogens in Pork 2001. Leipzig, Germany September 2 - 5, 2001
8. Jacquet C et al. 1995. Investigations related to the epidemic strain involved in the French listeriosis outbreak in 1992. *Appl. Environ. Microbiol.* 61: 2242 - 2246.
9. Murray, Patrick R., Ken S. Rosenthal, George S. Kobayashi, and Michael A. Pfaller. *Medical Microbiology*. 1998. St. Louis, MO.
10. Wesley IV. 1998. Listeriosis in animals. In *Listeria, Listeriosis and Food Safety*, ET Ryser and EH Marth (editors). pp. 39-74. Marcel Dekker, Inc., New York.

Table 1. Sequences and concentrations of primers and probe for amplification of the *iap* gene of *Listeria* species.

| Primer/Probe | Annealing site | Sequence | Concentration (nM) |
|-----------------|----------------|------------------------------------------|--------------------|
| <i>iap</i> -P-5 | 49 to 68 | 5' CGG TAA CAG CAT TTG CTG CGC CAA C 3' | 250 |
| <i>iap</i> -F | 1 to 26 | 5' ATG AAT ATG AAA AAA GCA ACT ATC GC 3' | 300 |
| <i>iap</i> -R | 1419 to 1437 | 5' TTA TAC GCG ACC GAA GCC A 3' | 50 |

Table 2. Specificity of *iap* assay for various *Listeria* strains.

| Bacteria | NADC strain | Amplification of <i>iap</i> |
|--------------------------------|-------------|-----------------------------|
| <i>L. monocytogenes</i> (4b) | 2045 | + |
| <i>L. monocytogenes</i> (1/2a) | 2052 | + |
| <i>L. monocytogenes</i> (4b) | 2840 | + |
| <i>L. innocua</i> | 2841 | + |
| <i>L. innocua</i> | 2842 | + |
| <i>L. murrayi</i> | 2843 | - |
| <i>L. seelegeri</i> | 2844 | + |
| <i>L. welshimeri</i> | 2845 | + |
| <i>L. ivanovii</i> | 6185 | + |
| <i>S. aureus</i> | 3649 | - |
| <i>S. aureus</i> | 3690 | - |
| <i>E. coli</i> | 8176 | - |
| <i>H. alvei</i> | 8526 | - |
| <i>P. vulgaris</i> | | - |
| <i>K. pneumoniae</i> | | - |
| <i>P. aeruginosa</i> | | - |
| <i>E. clocae</i> | | - |
| <i>C. coli</i> | 5095 | - |
| <i>E. coli</i> | 8252 | - |